

G-Force™ SureFit™



Taylor™

Healthcare Products, Inc.

DISPOSABLE FITTED STRETCHER SHEET SERIES

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's durable SureFit™ fitted stretcher sheets are the best disposable fitted sheets available! The Taylor G-Force™ Fitted Sheet Series is designed specifically for the Stryker® XPS mattress. Designed and tested to aid in transfer of patients up to 350 pounds, the Taylor G-Force™ Fitted Sheet Series has you covered.

- Custom fitted to the unique XPS mattress design
- Available in fluid resistant* fluid impervious** and heavy duty versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets
- The G-Force fitted sheets series is designed to aid in patient transfers up to 350 pounds



| Item # | Product | Color | Material | Size | Packaging |
|-------------|---|-------|------------------|-----------|-----------|
| 90-GFRC3690 | G-Force™ Fitted Sheet, for Stryker® XPS | Blue | Fluid Resistant | 36" x 90" | 30/case |
| 90-IFRC3690 | I-Force™ Fitted Sheet, for Stryker® XPS | Blue | Fluid Impervious | 36" x 90" | 30/case |
| 90-SUPER-G | SUPER-G™ Fitted Sheet, for Stryker® XPS | White | Fluid Resistant | 36" x 90" | 30/case |

*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure

**Fluid Impervious - Fluid does not penetrate the material even with time and pressure

LP90GFRCAD0417