



TITAN™ & TITAN II™ SOFT STRETCHERS

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service! The Titan Series are products no service can do without. They are essential to keep patients and rescuers safe as they navigate challenges faced daily.

- Cost effective, reusable, easy to store and durable
- The essential tool to aid in lifting bariatric patients
- Vital to have on the scene when space limitations prohibit use of standard EMS stretchers
- 16 SureGrip™ rubber handles ensure rescuers the comfort and stability of a firm grip, which provides easy lifting while transporting a patient
- Made of ultra-strong PVC and nylon webbing
- **TITAN** is designed to accommodate 16" – 18" backboards
- **TITAN II** is designed to be used without a backboard
- Easily disinfected and air dried
- Folds for space-saving storage in the accompanying custom Titan carrying case
- Available individually or in cases of 5
- Lifting capacity of 1500 lbs



Item #	Product	Color	Material	Size	Packaging
44-TT4080	Titan	Black	PVC, nylon webbing, SureGrip™ rubber handles	40"x 80"	5 per case
44-TITAN II	Titan II	Black	PVC, nylon webbing, SureGrip™ rubber handles	40"x 80"	5 per case

*Backboard not included

LP44A1018


Taylor®
Healthcare Products, Inc.