




# Yellow Emergency Blankets

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's Yellow Emergency Blankets are extremely durable and ideal to use when covering patients and/or equipment during recoveries and transports. The highly visible yellow color adds a layer of safety at accident scenes.

- Water proof
- Nonwoven PP material laminated to PE material which offers a heavy duty impervious barrier to the elements
- Retains body heat
- Color: Highly visible Yellow/White
- Available in two sizes: 60" x 90" and 50" x 84"



Item #	Product	Color	Material	Size	Packaging
90-YEB6090	Yellow Emergency Blanket	Yellow/White	Poly/Nonwoven	60" x 90"	20 per case
90-YEB5084	Yellow Emergency Blanket	Yellow/White	Poly/Nonwoven	50" x 84"	20 per case

LP90YEB0417