



TITAN^{XL} Soft Stretcher

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Our newest member of the Taylor Titan™ line of soft stretchers is here. Welcome the TitanXL™ to your bariatric unit, MCI protocol and disaster relief supply list. 20" wider than the standard patient carrier and 6" longer with the same exceptional Taylor quality. The extra width assists EMS professionals in keeping their hands clear of patient contact.

- Limited use soft stretchers
- The essential tool for lifting bariatric patients
- Vital to have when space limitations prohibit use of standard EMS stretchers
- 16 SureGrip™ rubber handles ensure rescuers the comfort & stability of a firm grip, providing easy lifting while transporting a patient
- Folds for space-saving storage in the accompanying custom TitanXL™ carrying case
- Lifting capacity of 1500 lbs.
- Fluid impervious* to help keep equipment clean & contain fluids
- Durable, compact, lightweight & cost effective
- Available individually or in cases of 5



Item #	Product	Color	Material	Size	Packaging
43-TITANXL	TitanXL™ Soft Stretcher	White	Nonwoven/Poly, SureGrip™ rubber handles	60" x 86"	5/case

*Fluid Impervious - Fluid does not penetrate the material even with time and pressure



LP43TXL102424